



# 2019-nCoV IgG & IgM

## A Diagnostic Solution for COVID-19

COVID-19 is a disease caused by the novel coronavirus (2019-nCoV). The disease can be screened and diagnosed by measuring the antibody IgG & IgM for this virus. Antibody IgM appears in the early stage of infection, while Antibody IgG in the late stage.


### Test kits

Novel Coronavirus (2019-nCoV) IgG/IgM Test Kit



### Features

- Colloidal Gold
- Qualitative testing with minimal interference
- 25T with individual pouch for each test
- Whole blood, serum and plasma
- Low sample volume (W.B. 20uL, S/P 10uL)
- Rapid results within 15 minutes
- Storage at room temperature



\* Take 2019-nCoV IgG/IgM test kit for example

### Interpretation

